

Factor IXa Control

Art No 9588

Lot# HKG2003

Reconstitute in 2.0 mL purified water

For In Vitro Laboratory Use Only – Not for Diagnostic Use

PROPERTIES

Lyophilized preparation of human Factor IXa. The Factor IXa Control is designed for direct use in the Rox FIX-A kit. No further dilution is necessary after reconstitution.

Factor IXa activity after reconstitution in **2.0 mL** purified water:

0.39 mIU/mL

Range (n=29): **0.34 – 0.45 mIU/mL**

The above Factor IXa activity apply only to **lot# HKF2003**.

The calibration has been made against the 1st IS FIXa 97/562 using the Rox FIX-A kit (Art No 95 00 30)

PACKAGE

10 x 2.0 mL

STORAGE

Keep at -20°C or below.

STABILITY

See expiry date on label

A reconstituted vial is stable for 8 hours at 2-8°C.

Rossix
Taljegårdsgatan 3B
S-431 53 Mölndal, Sweden

www.rossix.com
info@rossix.com