

Factor IXa Calibrator

Art No 9599

Lot# GSF2239

Reconstitute in 2.0 mL purified water

For Laboratory Use Only – Not for Diagnostic Use

PROPERTIES

Lyophilized preparation of human Factor IXa. Factor IXa activity after reconstitution in **2.0 mL** purified water:

6.1 mIU/mL

Range (95% CI): **5.9 – 6.4 mIU/mL**

The above Factor IXa activity apply only to **lot# GSF2239**.

The calibration has been made against the 1st IS FIXa 97/562 using the Rox FIX-A kit (Art No 95 00 30)

PACKAGE

10 x 2.0 mL

STORAGE

Keep at -20°C or below.

STABILITY

See expiry date on label

A reconstituted vial is stable for 8 hours at 2-8°C.