

Factor IXa Control

Art No 9588

Lot# HKF2239

Reconstitute in 2.0 mL purified water

For Laboratory Use Only – Not for Diagnostic Use

PROPERTIES

Lyophilized preparation of human Factor IXa. The Factor IXa Control is designed for direct use in the Rox FIX-A kit. No further dilution is necessary after reconstitution.

Factor IXa activity after reconstitution in **2.0 mL** purified water:

0.79 mIU/mL

Range (95% CI): **0.76 – 0.82 mIU/mL**

The above Factor IXa activity apply only to **lot# HKF2239**.

The calibration has been made against the 1st IS FIXa 97/562 using the Rox FIX-A kit (Art No 95 00 20)

PACKAGE

10 x 2.0 mL

STORAGE

Keep at -20°C or below.

STABILITY

See expiry date on label

A reconstituted vial is stable for 8 hours at 2-8°C.

Rossix
Taljegårdsgatan 3B
S-431 53 Mölndal, Sweden

www.rossix.com
info@rossix.com